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[Menu] → [Setting] → [Assay Parameters] → [PT THS] → [Test Protocol]

Assay Group	Basic	Re-analysis	Reflex	QC	Test Protocol			
PT THS	Rinse at the Time of Primary							
APTT PSL	Aspiration	<input type="checkbox"/>	Pre-Rinse	<input type="checkbox"/>	Post-Rinse			
Fibrinogen	Sample Dispensing /Dilution							
Thrombine Time	Dilution Setting							
AT III	Buffer:	OVB	Dilution Ratio	1/1, 1/2, 1/4, 1/8		Edit	Add	
BCPC	Aspirated Sample	Sample Aspiration Vol [μl]	Buffer Volume [μl]	Pre-rinse List	Post-rinse List	Edit		
F VIII	Sample	50	0	-	-	Delete		
F IX								
BCPIg								
	Factor-Deficient Plasma							
	Aspirated Reagent	Reagent Aspiration Vol [μl]	Pre-rinse List	Post-rinse List	Edit			
	Add							
	Aspirated Reagent	Reagent Aspiration Vol [μl]	Pipette at [sec]	Mixing Level	Pre-rinse List	Post-Rinse List	Edit	
	PT THS	100	190	Weak	-	Clean I	Delete	
	Measurement Condition							
	Measurement Time (Main)	180 sec	Measurement Time (Sub)	180 sec	Seq.	Reagent	Lot management	Stable
					1	OVB	<input type="checkbox"/>	0
					2	PT THS	<input checked="" type="checkbox"/>	0
	Measurement Time (Main)	Measurement Time (Sub)						
	<input type="checkbox"/>	Mixing during measurement						
	STEP 1	0 sec –	0	0 rpm				
	STEP 2	0 sec –	0	0 rpm				
	STEP 3	0 sec –	180 sec	0 rpm				
	Export	Add	Load	Print	Save	Close		

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[Menu] → [Setting] → [Assay Parameters] → [Fbg] → [Test Protocol]

Assay Group	Basic	Re-analysis	Reflex	QC	Test Protocol
PT THS	Rinse at the Time of Primary				
APTT PSL	Aspiration <input type="checkbox"/> Pre-Rinse <input type="checkbox"/> Post-Rinse				
Fibrinogen	Sample Dispensing /Dilution				
Thrombine Time	Dilution Setting				
AT III	Buffer:	OVB	Dilution Ratio	1/1, 1/2, 1/4, 1/8, 2/1, 3/1	Edit Add
BCPC	Aspirated Sample	Sample Aspiration Vol [μl]	Buffer Volume [μl]	Pre-rinse List	Post-rinse List Edit
FVIII	Sample	10	90	-	- Delete
FIX					
BCPIg					
Factor-Deficient Plasma					Add
Aspirated Reagent					Reagent Aspiration Vol [μl]
					Pre-rinse List Post-rinse List Edit
					Delete
					Add
Aspirated Reagent					Reagent Aspiration Vol [μl]
Fbg					50 Pipette at [sec] 190 Mixing Level Weak Pre-rinse List - Post-Rinse List Clean I Edit
					Delete
Measurement Condition					
Measurement Time (Main)		100 sec	Measurement Time (Sub)	100 sec	Seq. Reagent Lot management Stable
					1 OVB <input type="checkbox"/> 0
Measurement Time (Main)		Measurement Time (Sub)			2 Fbg <input checked="" type="checkbox"/> 0
<input type="checkbox"/>		Mixing during measurement			
STEP 1	0 sec –	0	0 rpm		
STEP 2	0 sec –	0	0 rpm		
STEP 3	0 sec –	100 sec	0 rpm		
Export					Add Load Print Save Close

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[Menu] → [Setting] → [Assay Parameters] → [ТТ] → [Test Protocol]

Assay Group	Basic	Re-analysis	Reflex	QC	Test Protocol									
PT THS	Rinse at the Time of Primary													
APTT PSL	Aspiration <input type="checkbox"/> Pre-Rinse <input type="checkbox"/> Post-Rinse													
Fibrinogen	Sample Dispensing /Dilution													
Thrombine Time	Dilution Setting													
AT III	Buffer:	OVB	Dilution Ratio	1/1	Edit	Add								
BCPC	Aspirated Sample	Sample Aspiration Vol [μl]	Buffer Volume [μl]	Pre-rinse List	Post-rinse List	Edit								
FVIII	Sample	50	0	-	-	Delete								
FIX														
BCPIg														
	Factor-Deficient Plasma						Add							
	Aspirated Reagent	Reagent Aspiration Vol [μl]	Pre-rinse List	Post-rinse List		Edit								
						Delete								
						Add								
	Aspirated Reagent	Reagent Aspiration Vol [μl]	Pipette at [sec]	Mixing Level	Pre-rinse List	Post-Rinse List	Edit							
	Test Thr	100	70	Weak	-	Clean I	Delete							
	Measurement Condition													
	Measurement Time (Main)	160 sec	Measurement Time (Sub)	160 sec	Seq.	Reagent	Lot management	Stable						
					1	OVB	<input type="checkbox"/>	0						
					2	Test Thr	<input checked="" type="checkbox"/>	0						
	Measurement Time (Main)	Measurement Time (Sub)												
	<input type="checkbox"/>	Mixing during measurement												
	STEP 1	0 sec –	0	0 rpm										
	STEP 2	0 sec –	0	0 rpm										
	STEP 3	0 sec –	160 sec	0 rpm										
	<table border="1"> <tr> <td>Export</td> <td>Add</td> <td>Load</td> <td>Print</td> <td>Save</td> <td>Close</td> </tr> </table>								Export	Add	Load	Print	Save	Close
Export	Add	Load	Print	Save	Close									

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[Menu] → [Setting] → [Assay Parameters] → [AT3] → [Test Protocol]

Assay Group	Basic	Re-analysis	Reflex	QC	Test Protocol
PT THS	Rinse at the Time of Primary				
APTT PSL	Aspiration <input type="checkbox"/> Pre-Rinse <input type="checkbox"/> Post-Rinse				
Fibrinogen	Sample Dispensing /Dilution				
Thrombine Time	Dilution Setting				
AT III	Buffer:	OVB	Dilution Ratio	1/1, 1/2, 1/4, 1/8, 3/2, 0/1	
BCPC	Aspirated Sample	Sample Aspiration Vol [μl]	Buffer Volume [μl]	Pre-rinse List	Post-rinse List
FVIII	Sample	21	147	-	-
FIX	Dilution	24	0		
BCPIg	Factor-Deficient Plasma				
	Aspirated Reagent	Reagent Aspiration Vol [μl]	Pre-rinse List	Post-rinse List	
	Aspirated Reagent	Reagent Aspiration Vol [μl]	Pipette at [sec]	Mixing Level	Pre-rinse List
	AT3thro	175	40	Weak	Clean I
	AT3Subs	33	230	Strong	
	Measurement Condition				
	Measurement Time (Main)	60 sec	Measurement Time (Sub)	60 sec	Seq.
					Reagent
					Lot management
					Stable
	Measurement Time (Main)	Measurement Time (Sub)			
	<input type="checkbox"/>	Mixing during measurement			
	STEP 1	0 sec – 0	0 rpm		
	STEP 2	0 sec – 0	0 rpm		
	STEP 3	0 sec – 60 sec	0 rpm		
	Export	Add	Load	Print	Save
					Close

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[Menu] → [Setting] → [Assay Parameters] → [BCPIg] → [Test Protocol]

Assay Group	Basic	Re-analysis	Reflex	QC	Test Protocol	
PT THS	Rinse at the Time of Primary					
APTT PSL	Aspiration <input type="checkbox"/> Pre-Rinse <input type="checkbox"/> Post-Rinse					
Fibrinogen	Sample Dispensing /Dilution					
Thrombine Time	Dilution Setting					
AT III	Buffer:	OVB	Dilution Ratio	1/1, 1/2, 1/4, 1/8, 3/2, 0/1		Edit Add
BCPC	Aspirated Sample	Sample Aspiration Vol [μl]	Buffer Volume [μl]	Pre-rinse List	Post-rinse List	Edit
FVIII	Sample	21	147	-	-	Delete
FIX	Dilution	28	0			
BCPIg	Factor-Deficient Plasma					
	Aspirated Reagent	Reagent Aspiration Vol [μl]	Pre-rinse List	Post-rinse List		Edit
	PIReag					Delete
	Add					
	Aspirated Reagent	Reagent Aspiration Vol [μl]	Pipette at [sec]	Mixing Level	Pre-rinse List	Post-Rinse List Edit
	PIReag	175	40	Middle	Clean I	Clean I Delete
	PISubs	35	350	Strong		
	Measurement Condition					
	Measurement Time (Main)	120 sec	Measurement Time (Sub)	120 sec	Seq.	Reagent Lot management Stable
					1	OVB <input type="checkbox"/> 0
					2	PIReag <input checked="" type="checkbox"/> 0
					3	PISubs <input checked="" type="checkbox"/> 0
	Measurement Time (Main)	Measurement Time (Sub)				
	<input type="checkbox"/>	Mixing during measurement				
	STEP 1	0 sec –	0	0 rpm		
	STEP 2	0 sec –	0	0 rpm		
	STEP 3	0 sec –	120 sec	0 rpm		
	Export Add Load Print Save Close					

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[Menu] → [Setting] → [Assay Parameters] → [BCPC] → [Test Protocol]

Assay Group	Basic	Re-analysis	Reflex	QC	Test Protocol			
PT THS	Rinse at the Time of Primary							
APTT PSL	Aspiration	<input type="checkbox"/>	Pre-Rinse	<input type="checkbox"/>	Post-Rinse			
Fibrinogen	Sample Dispensing							
Thrombine Time	/Dilution							
AT III	Dilution Setting							
	Buffer:	OVB	Dilution Ratio	1/1, 1/2, 1/4, 1/8, 3/2, 0/1			Edit Add	
BCPC	Aspirated Sample	Sample Aspiration Vol [μl]	Buffer Volume [μl]	Pre-rinse List	Post-rinse List	Edit		
FVIII	Sample	15	15	-	-	Delete		
FIX								
BCPIg								
	Factor-Deficient Plasma							
	Aspirated Reagent	Reagent Aspiration Vol [μl]	Pre-rinse List	Post-rinse List	Edit			
					Delete			
	Aspirated Reagent	Reagent Aspiration Vol [μl]	Pipette at [sec]	Mixing Level	Pre-rinse List	Post-Rinse List	Edit	
	BCPCAct	150	60	Weak	Clean I	Clean I	Delete	
	BCPCSubs	30	350	Strong				
	Measurement Condition							
	Measurement Time (Main)	80 sec	Measurement Time (Sub)	80 sec	Seq.	Reagent	Lot management	Stable
					1	OVB	<input type="checkbox"/>	0
	Measurement Time (Main)	Measurement Time (Sub)			2	BCPCAct	<input checked="" type="checkbox"/>	0
	<input type="checkbox"/>	Mixing during measurement			3	BCPCSubs	<input checked="" type="checkbox"/>	0
	STEP 1	0 sec –	0	0 rpm				
	STEP 2	0 sec –	0	0 rpm				
	STEP 3	0 sec –	80 sec	0 rpm				
	Export	Add	Load	Print	Save	Close		

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[Menu] → [Setting] → [Assay Parameters] → [LA R] → [Test Protocol]

BA-TEST

Assay Group	Basic	Re-analysis	Reflex	QC	Test Protocol			
PT THS	Rinse at the Time of Primary							
APTT PSL	Aspiration	<input type="checkbox"/>	Pre-Rinse	<input type="checkbox"/>	Post-Rinse			
Fibrinogen	Sample Dispensing /Dilution							
Thrombine Time	Dilution Setting							
AT III	Buffer:	OVB	Dilution Ratio	1/1, 1/2, 1/4, 1/8, 1/16, 1/32, 3/2, 2/1, 1/64, 1/256			Edit Add	
BCPC	Aspirated Sample	Sample Aspiration Vol [μl]	Buffer Volume [μl]	Pre-rinse List	Post-rinse List	Edit		
LA-1 R	Sample	50	-	-	-	Delete		
LA 2 R	Dilution	-	-					
BCPIg								
	Factor-Deficient Plasma							
	Aspirated Reagent	Reagent Aspiration Vol [μl]	Pre-rinse List	Post-rinse List			Edit	
							Delete	
							Add	
	Aspirated Reagent	Reagent Aspiration Vol [μl]	Pipette at [sec]	Mixing Level	Pre-rinse List	Post-Rinse List	Edit	
	LA-1	50	20	Weak			Delete	
	CaCl2	50	90	Strong				
	Measurement Condition							
	Measurement Time (Main)	360 sec	Measurement Time (Sub)	360 sec	Seq.	Reagent	Lot management	Stable
	<input type="checkbox"/>	Mixing during measurement			1	OVB	<input type="checkbox"/>	0
	STEP 1	0 sec –	0	0 rpm	3	LA-1	<input checked="" type="checkbox"/>	0
	STEP 2	0 sec –	0	0 rpm	4	CaCl2	<input type="checkbox"/>	
	STEP 3	0 sec –	360 sec	0 rpm				
	Export	Add	Load	Print	Save	Close		

Assay Group	Basic	Re-analysis	Reflex	QC	Test Protocol			
PT THS	Rinse at the Time of Primary							
APTT PSL	Aspiration <input type="checkbox"/> Pre-Rinse <input type="checkbox"/> Post-Rinse							
Fibrinogen	Sample							
Thrombine Time	Dispensing /Dilution							
AT III	Dilution Setting	OVB	Dilution Ratio	1/1, 1/2, 1/4, 1/8, 1/16, 1/32, 3/2, 2/1, 1/64, 1/256			Edit Add	
BCPC	Aspirated Sample	Sample Aspiration Vol [μl]	Buffer Volume [μl]	Pre-rinse List	Post-rinse List		Edit	
LA-1 R	Sample	50	-	-	-		Delete	
LA 2 R	Dilution	-	-					
BCPIg								
	Factor-Deficient Plasma							Add
	Aspirated Reagent	Reagent Aspiration Vol [μl]	Pre-rinse List	Post-rinse List			Edit	
							Delete	
								Add
	Aspirated Reagent	Reagent Aspiration Vol [μl]	Pipette at [sec]	Mixing Level	Pre-rinse List	Post-Rinse List	Edit	
	LA-2	50	20	Weak			Delete	
	CaCl2	50	90	Strong				
	Measurement Condition							
	Measurement Time (Main)	360 sec	Measurement Time (Sub)	360 sec	Seq.	Reagent	Lot management	Stable
	<input type="checkbox"/>	Mixing during measurement			1	OVB	<input type="checkbox"/>	0
					3	LA-2	<input checked="" type="checkbox"/>	0
					4	CaCl2	<input type="checkbox"/>	
	STEP 1	0 sec	0	0 rpm				
	STEP 2	0 sec	0	0 rpm				
	STEP 3	0 sec	360 sec	0 rpm				
	Export	Add	Load	Print	Save	Close		

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**Assay Group Setup**

Assay Group: **Red D-dimer**

Basic | Re-analysis | Reflex | QC | Test Protocol

Rinse at the Time of Primary Aspiration:  Pre-Rinse  Post-Rinse

Sample Dispensing/Dilution

Dilution Setting: Buffer: **DD Diluent** Dilution Ratio for Measurements: 1/1, 1/2, 1/4, 1/8, 1/16, 1/32

Aspirated Sample	Sample Aspiration Vol. [uL]	Buffer Volume [uL]	Pre-Rinse List	Post Rinse List
Sample	10	10	-	Clean 1

Factor-Deficient Plasma

Aspirated Reagent	Reagent Aspiration Vol.[uL]	Pre-Rinse List	Post Rinse List

Reagent Protocol

Aspirated Reagent	Reagent Aspiration Vol.[uL]	Pipette at [sec]	Mixing Level	Pre-Rinse List	Post-Rinse List
DD RBuff	80	20	Middle	-	Clean 1
DD Latex	80	90	Middle	-	Clean 1

Measurement Condition

Measurement Time(Main): 200 sec Measurement Time (Sub): 200 sec

Mixing during measurement:

STEP1: 0 sec - 0 rpm

STEP2: 0 sec - 0 rpm

STEP3: 0 sec - 200 sec 0 rpm

Seq.	Reagent	Lot Manag	Stable Time
1	DD Diluent	<input type="checkbox"/>	0
2	DD RBuff	<input type="checkbox"/>	0
3	DD Latex	<input type="checkbox"/>	0

Buttons: Copy, Up, Delete, Down, Export, Add, Load, Print, Save, Close

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Graph Axis  
X Axis Min.:     X Axis Max.:     Y Axis Min.:     Y Axis Max.:

Dilution Analysis:

Correction	Calibrator	Dilution Ra	Replicatio	Stable Time
------------	------------	-------------	------------	-------------

Assay Parameter Setup

Assay Parameter:       Display Name:       Management ID:

Valid      Updated on:       Assay Parameter Host ID:

Calculation Method:  Evaluation Preset

Report Limit Check

Upper Limit:       Lower Limit:

Mark Limit Check

Upper Limit:       Lower Limit:

Replication Difference Limit Check

Upper Limit:

MDA Slope Ratio Check

Upper Limit:       Lower Limit:

Assay Parameter Setup

Assay Parameter:  Display Name:  Management ID:

Valid Updated on:  Assay Parameter Host ID:

Calculation Method | Data Check | Evaluation Preset

Evaluation Algorithm:  Detection Principle:

Target

Wave:

Type:

Gain:

Smoothing

Reference Correction

Median Preceding Point:  Following Point:

Moving Average Preceding Point:  Following Point:

Correction

User

A (Slope)	B (Offset)
<input type="text" value="1.00"/>	<input type="text" value="0.0000"/>

Evaluation Parameter | Evaluation Check Parameter

Start Time[sec]:  Polynomial Dimension:

End Time[sec]:  Integral Area(auc):

Min. Regression Time[sec]:

Variable Start Time

Preeval Search Width[sec]:	<input type="text" value="50"/>
Upper Eval Offset[sec]:	<input type="text" value="60"/>
Upper Preeval Rate[dOD/min]:	<input type="text" value="0.0000"/>
Lower Preeval Rate[dOD/min]:	<input type="text" value="0.1000"/>
Minimum Search Window[sec]:	<input type="text" value="20"/>

Save Close

Assay Parameter Setup

Assay Parameter:  Display Name:  Management ID:

Valid Updated on:  Assay Parameter Host ID:

Calculation Method | Data Check | Evaluation Preset

Evaluation Algorithm:  Detection Principle:

Target  
Wave:   
Type:   
Gain:

Smoothing  
 Reference Correction  
 Median Preceding Point:  Following Point:   
 Moving Average Preceding Point:  Following Point:

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Mark Limit Check

Upper Limit:       Lower Limit:

Replication Difference Limit Check

Upper Limit:

MDA Slope Ratio Check

Upper Limit:       Lower Limit:

Save      Close

Assay Parameter Setup

Assay Parameter:  Display Name:  Management ID:

Valid Updated on:  Assay Parameter Host ID:

Calculation Method | Data Check | Evaluation Preset

Evaluation Algorithm:  Detection Principle:

Target

Wave:

Type:

Gain:

Smoothing

Reference Correction

Median Preceding Point:  Following Point:

Moving Average Preceding Point:  Following Point:

Correction

User

A (Slope)	B(Offset)
<input type="text" value="1.00"/>	<input type="text" value="0.0000"/>

Evaluation Parameter | Evaluation Check Parameter

Start Time[sec]:  Polynomial Dimension:

End Time[sec]:  Integral Area(auc):

Min. Regression Time[sec]:

Variable Start Time

Preeval Search Width[sec]:	<input type="text" value="50"/>
Upper Eval Offset[sec]:	<input type="text" value="60"/>
Upper Preeval Rate[dOD/min]:	<input type="text" value="0.0000"/>
Lower Preeval Rate[dOD/min]:	<input type="text" value="0.1000"/>
Minimum Search Window[sec]:	<input type="text" value="20"/>

Save Close

Assay Parameter Setup

Assay Parameter:  Display Name:  Management ID:

Valid Updated on:  Assay Parameter Host ID:

Calculation Method | Data Check | Evaluation Preset

Evaluation Algorithm:  Detection Principle:

Target  
Wave:   
Type:   
Gain:

Smoothing  
 Reference Correction  
 Median Preceding Point:  Following Point:   
 Moving Average Preceding Point:  Following Point:

Correction  
 User  
A (Slope) B (Offset)

Evaluation Parameter | Evaluation Check Parameter

High Light Limit:   
Low Light Limit:

Reaction Curve  
Slope:   
Slope Thresh.:

Antigen Check  
Start Time 1 [sec]:   
End Time 1 [sec]:   
Cut Off:   
Rate:   
Offset:   
Start Time 2 [sec]:   
End Time 2 [sec]:

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